

Clean Room Solutions-Biological Safety and Automations

Biosafety Cabinet

- Laminar Air Flow
- Sampling Booth
- PCR Workstation
- Scrub Station



Technology To Enable



Biosafety Controller with Differential Sensor



L: 310mm x W: 95mm x H: 50mm

- Model:**
- Dual Sensor (Duo)
 - Single Sensor (Uni)
 - Without Sensor

- Password (Programmable) for Authorized Operation.
- Built in differential Pressure display (mmWC / Pascal).
- Air Velocity Display and set the safety parameters.
- Initial Self Cleaning Cycle.
- Real Time Clock.
- UV Timer.
- Total UV used Time.
- Total HEPA Filter run time.
- HEPA Filter Alarm after programmed hours with ON/OFF. Etc.




Biosafety Controller with Air Velocity and Differential Pressure



L: 310mm x W: 95mm x H: 50mm

- Model:**
- Dual Sensor (Duo)
 - Single Sensor (Uni)
 - Without Sensor

- Password (Programmable) for Authorized Operation.
- Built in differential Pressure display (mmWC / Pascal).
- Air Velocity Display and set the safety parameters.
- Initial Self Cleaning Cycle.
- Real Time Clock.
- UV Timer.
- Total UV used Time.
- Total HEPA Filter run time.
- HEPA Filter Alarm after programmed hours with ON/OFF. Etc.



Laminar Air Flow Controller with Differential Sensor Product Deck



L: 310mm x W: 95mm x H: 50mm

- Model:**
- Dual Sensor (Duo)
 - Single Sensor (Uni)
 - Without Sensor

- Password (Programmable) for Authorized Operation.
- Built in differential Pressure display (mmWC / Pascal).
- Air Velocity Display and set the safety parameters.
- Initial Self Cleaning Cycle.
- Real Time Clock.
- UV Timer.
- Total UV used Time.
- Total HEPA Filter run time.
- HEPA Filter Alarm after programmed hours with ON/OFF. Etc.

